



LIMING BIO

®

Strongstep®

SARS-CoV-2 Antigen Rapid Test

A Community of Shared Future for Mankind

Introduction of Company

Nanjing Liming Bio-Products Co., Ltd

It is a professional manufacturer of in vitro diagnostic reagents

2001

Time of establishment

10%

Percentage of master degree or above

106

Existing employees

1500

Purification plant




Strongstep®

- » Rapid Diagnostic Test
- » 5 million Tests per Year

Limingbio diagnostics are ASSURED


**SIMPLE
SAFE
FAST
ACCURATE**



Strongstep® SARS-CoV-2 Antigen Rapid Test Cassette



22.5×18×6.5CM



REF:500200

20 Tests/Kit



Specimens: nasal swab / nasopharyngeal swab

20 Individually packed test devices

2 Extraction Buffer vials

20 Extraction tubes


1 Workstation

1 Package insert

What is the Rapid Antigen Test ?

A novel coronavirus antigen was detected by latex immunochromatography in the nasal swabs and nasopharyngeal swabs.

Sampling illustration




nasal swab
nasopharyngeal swab

INSTRUCTION to USE


Sample Treatment Solution



Squeeze the tube and rotate the swab 15 times



Remove the test card from the foil bag
Add 3 drops of specimen solution to round sampling hole



Waiting for 15 minutes
reading result

Result



positive



negative



invalid



Strongstep® SARS-CoV-2 Antigen Rapid Test System Device



REF:500210

20 Tests/Kit




Specimens: nasal swab / nasopharyngeal swab

CONTENTS


- ◆ Individually Packed Test Devices : 20
Each device contains a strip with colored conjugates and reactive reagents pre-spread at the corresponding regions.
- ◆ Extraction Buffer vials in the device : 20
0.1 M Phosphate buffered saline (PBS) and 0.02% sodium azide.
- ◆ Package insert : 1
- ◆ Workstation : 1



INSTRUCTION to USE



Invert a sterile swab into the nostril of the patient, reaching the surface of the posterior nasopharynx. Swab over the surface of the posterior nasopharynx.



withdraw the sterile swab from the nasal cavity




Unscrew the cover of the device



×5




Put the swab into the tube, break the swab with the breakpoint, let the sampled swab fall into the tube and discard the upper stick




Screw the cover of the device




Break the blue stick.



FIRMLY squeeze the blue tube, make sure all the liquid fall into the lower tube




Vortex the device vigorously



Invert the device, let the sample buffer migrate onto the test strip




Put the device into the Plastic base



Plastic base




Read the results at 15 minutes



Advantages of midstream:

- 1、 Double biosafety protection design to protect operator and Lab
- 2、 Independent Packaging
- 3、 Can be sold in supermarket or drugstores, can be used for home self-testing

INTERPRETATION OF RESULTS

POSITIVE RESULT: 	Two colored bands appear within 15minutes. One colored band appears in the Control Zone (C) and another colored band appears in the Test Zone (T). The test result is positive and valid. No matter how faint the colored band appears in the Test Zone (T), the test result should be considered a s positive result.
NEGATIVE RESULT: 	One colored bands appears in the Control Zone (C) within 15 minutes. No colored band appears in the Test Zone (T). The test result is negative and valid.
INVALID RESULT: 	No colored band appears in the Control Zone (C) within 15 minutes. The test result is invalid. Repeat the test with a new test device.

Performance

Positive Percent Agreement:


(PPA)=73.02% (60.35% ~ 83.43%)

Negative Percent Agreement:

(NPA)=98.08% (93.23% ~ 99.77%)

	PCR Comparator			
StrongStep [®] SARS CoV-2 Antigen Test	Antigen Test	Positive	Negative	Total
	Positive	46	2	48
	Negative	17	102	119
	Total	63	104	167

LIMING issues the "List Of COVID-19 In Vitro Diagnostic Manufacturers and Exporters" and updates it regularly on the web site.



ISO 13485



Medicines & Healthcare products
Regulatory Agency



Please use RG3, the Registration form, to tell us about any of these changes. A fee of £70 is payable for each change or set of changes notified.

Thank you for registering the following generic groups of devices

1. **Part 5: IVDs which are not Annex II and not self-test devices**
- 2.
3. **For reagents, reagent products, calibration and control materials:**
4. **group by common technological characteristics and/or analytes**
- 5.
6. **New products:**
7. **None**
- 8.
9. **For performance evaluation:**
10. **None**
- 11.
12. **Neither:**
13. **HSV Antigen**
14. **Gonococcal Antigen Detection**
15. **Candida albicans**
16. **Other Parasitology**
17. **Other Bacteriology Rapid Tests**
18. **Rotavirus**
19. **Adenovirus**
20. **Other Multiple Viruses**
21. **H. Pylori Antibody Assays**
22. **H. Pylori Antigen Detection**
23. **Strep B - Rapid Test**
24. **Human Papilloma Virus**
25. **Strep A - Rapid Test**
26. **Haemoglobin (Hb)**
27. **Procalcitonin**
28. **Salmonella Antigen Detection**
29. **Salmonella Antibody Assays**
30. **Legionella Antibody Assays**
31. **Other Mycology Immunoassays**
32. **Other Specific Proteins Rapid Tests**
33. **Other Individual and Specified Hormones/Proteins RT & POC**
34. **Coronavirus**
35. **Coronavirus - NA Reagents**
- 36.
- 37.
38. **For other IVDs, group by appropriate indications**
- 39.
40. **New products**

(MHRA) CE certificate

COVID-19 outbreak diagnostic

https://www.finddx.org/covid-19/

FIND
Because diagnosis matters

WHO WE ARE WHAT WE DO RESOURCES NEWSROOM PARTNERS & DONORS

CONTACT US

- [PerfectQ COVID-19 Coronavirus Real Time PCR Kit \(RUO\)](#)
- [JN Medsys, ProTect Covid-19 RT-qPCR kit \(RUO\) Contact](#)
- [KH Medical Co., Ltd, RADI COVID-19 Detection Kit and RADI COVID-19 Triple Detection Kit \(RUO\) Contact](#)
- [KogeneBiotech Co., Ltd, PowerChekTM 2019-nCoV Real-time PCR Kit \(Korea MFDS-EUA; CE-IVD\) Contact](#)
- [Lifiver, Lifiver Novel Coronavirus \(2019-nCoV\) Real Time Multiplex RT-PCRT kit \(China FDA-EUA; CE-IVD\)](#)
- [Liming Bio-Products Co., Ltd, *SrongStep® Novel Coronavirus \(SARS-CoV-2\) Multiplex Real-Time PCR Kit \(CE-IVD\) Contact*](#)
- [Luminex Corp., NxTAG CoV Extended Panel \(RUO\) Contact](#)
- [Mabsky Bio-Tech Co., Ltd Real-Time PCR Method Contact](#)
- [COVID-19 virus \(2019-nCoV\) Dual-Detection Kit \(RUO\)](#)
- [Influenza A virus, Influenza B virus & COVID-19 virus \(2019-nCoV\) Triple-Detection Kit \(RUO\)](#)
- [Medical Innovation Ventures Sdn Bhd, GenoAmp® Real-Time RT-PCR SARS-CoV-2 \(RUO\) Contact](#)
- [Nanjing Yazyne Medical Technology Co., LTD., 2019-Novel Coronavirus \(2019-nCoV\) Triplex RT-qPCR Detection Kit \(CE-IVD\) Contact](#)
- [NanoBio Lab, A*STAR Research Entities, Isothermal Exponential Amplification for COVID-19 Detection \(RUO\) Contact](#)
- [National Institute for Control of Vaccines and Biologicals, Accupid nCoV 2019 Detection Kit \(RUO\) Contact](#)
- [Ningbo Health Gene Technologies Co., Ltd, SARS-CoV-2 Virus Detection Diagnostic Kit \(RT-qPCR\)](#)

COVID-19 outbreak diagnostic

https://www.finddx.org/covid-19/

FIND
Because diagnosis matters

WHO WE ARE WHAT WE DO RESOURCES NEWSROOM PARTNERS & DONORS

CONTACT US

- [Anti-SARS-CoV-2 ELISA \(IgG\) \(manual; automated; RUO\)](#)
- [Guangzhou Darui Biotechnology Co., Ltd Contact](#)
- [2019 Novel Coronavirus \(2019-nCoV\) IgM Antibody Detection Kit \(ELISA Method\) \(RUO\)](#)
- [2019 Novel Coronavirus \(2019-nCoV\) IgG Antibody Detection Kit \(ELISA Method\) \(RUO\)](#)
- [Novel Coronavirus 2019-nCoV IgM Antibody Detection Kit \(Colloidal Gold Method\) \(RUO\)](#)
- [Novel Coronavirus 2019-nCoV IgG Antibody Detection Kit \(Colloidal Gold Method\) \(RUO\)](#)
- [Guangzhou Wondfo Biotech Co., Ltd, Finecare SARS-CoV-2 Antibody Test \(manual; RUO\) Contact](#)
- [Liming Bio-Products Co., Ltd, COVID-19 Antigen Rapid Test Device \(CE-IVD\) Contact](#)
- [SD BIOSENSOR Inc., STANDARD F COVID-19 Ag IIA \(manual; CE-IVD\) Contact](#)
- [Shenzhen Yhlo Biotech Co., Ltd Contact](#)
- [iFlash-SARS-CoV-2 IgM \(CE-IVD\)](#)
- [iFlash-SARS-CoV-2 IgG \(CE-IVD\)](#)
- [Taizhou ZECEN Biotech Co., Ltd., Contact](#)
- [SARS-CoV-2 IgM \(CE-IVD\)](#)
- [SARS-CoV-2 IgG \(CE-IVD\)](#)
- [Sugentech, Inc. \(manual, RUO\) Contact](#)
- [SGTi-flex COVID-19 IgM/IgG](#)
- [SGTi-flex COVID-19 IgM](#)
- [SGTi-flex COVID-19 IgG](#)

Access to WHO
emergency approval channel

Daftar Rekomendasi RDT Antibodi COVID-19

Diperbaharui tanggal 28 April 2020
Gugus Tugas Percepatan Penanganan COVID-19

Rekomendasi Merk RDT Antibodi COVID-19 yang dapat digunakan di Indonesia, diklasifikasikan sebagai berikut:

- **Rekomendasi:** Berdasarkan edaran WHO (WHO *does not endorse any of these products*) yang dipublikasikan pada 27 April 2020 terkait alat diagnostik untuk Diagnosis COVID-19 dan Sertifikasi oleh CE (Sertifikasi yang dikeluarkan Uni Eropa), FDA (Sertifikasi yang dikeluarkan Amerika Serikat), atau sertifikasi yang setara.
- **Alternatif Rekomendasi:** Berdasarkan edaran WHO yang dipublikasikan pada 27 April 2020 terkait alat diagnostik untuk Diagnosis COVID-19 atau Sertifikasi yang disebut di atas


Daftar Rekomendasi RDT Antibodi COVID-19

Diperbaharui tanggal 28 April 2020
Gugus Tugas Percepatan Penanganan COVID-19

Alternatif Rekomendasi RDT Antibodi untuk COVID-19

31. Covid-19 IgG/IgM Antibody Rapid Test Kit - Wuhan UNscience Biotechnology Co., Ltd
32. COVID-19 IgG/IgM BIO - QUIBASA QUIMICA BASICA LTDA
33. COVID-19 IgG/IgM Combo Rapid Test Device - Liming Bio-Products Co., Ltd
34. COVID-19 IgG/IgM Detection Kit (Colloidal Gold) - Hunan Lituo Biotechnology Co., Ltd
35. COVID-19 IgG/IgM Duo (automated) - NanoEnTek
36. COVID-19 IgG/IgM ECO Test - Eco Diagnostica Ltd
37. COVID-19 IgG/IgM 15 - ADVAGEN BIOTECH LTDA

Access to Indonesian official
website procurement list



A large number of orders have been received from the United States, South Korea, Japan, Italy, Hungary, Austria, Saudi Arabia, Iran, France, Spain, Georgia, Chile, Brazil, the Philippines, Mexico, Bulgaria and other countries.

We Shall Overcome Together

A Community of Shared Future for Mankind



LIMING BIO

THANKS