



LIMING BIO
diagnostics are ASSURED

www.stddiagnostics.com www.limingbio.com

*Specialized In Manufacturing
Rapid Tests For
Sexually Transmitted Diseases*

StrongStep® Gonorrhoeae/Chlamydia Antigen Combo Rapid Test



An immunochromatographic assay for qualitative presumptive detection of *Neisseria gonorrhoeae*/*Chlamydia trachomatis* antigen from female cervical swab and male urethral swab

One test for detecting *Gonorrhoeae* and *Chlamydia* simultaneously

- ◆ **Accurate:** 97.51% sensitivity, 97.40% specificity for *Gonorrhoeae*, 94.39% sensitivity, 99.55% specificity for *Chlamydia* compared with Culture
- ◆ **Rapid:** one test for two organisms with the same specimen in just 15 minutes
- ◆ **Early detection:** stem the development of PID and other complications at early stage
- ◆ **User-friendly:** easily performed and interpreted by all health care personnel
- ◆ **Room temperature storage**

Cat.	Product
500020	StrongStep® <i>Neisseria Gonorrhoeae</i> Antigen Rapid Test
500010	Strong Step® <i>Chlamydia Trachomatis</i> Antigen Rapid Test
500050	Strong Step® <i>Gonorrhoeae/Chlamydia</i> Antigen Combo Rapid Test

Detection	Specimen
Antigen	Cervical/urethral swab
Antigen	Cervical/urethral swab
Antigen	Cervical/urethral swab