



LIMING BIO  
diagnostics are ASSURED

[www.stddiagnostics.com](http://www.stddiagnostics.com) [www.limingbio.com](http://www.limingbio.com)

*Specialized In Manufacturing  
Rapid Tests For  
Sexually Transmitted Diseases*

**A Breakthrough in the Virologic Diagnosis of HSV Infection**

## StrongStep® HSV 1/2 Antigen Rapid Test



**An immunochromatographic assay for qualitative presumptive detection of HSV 1 and 2 antigen from Mucocutaneous lesions**

- More accurate and specific than the type-specific HSV serologic assays (IgM/IgG) in symptomatic patients.
- Can be used in the differential diagnosis of genital ulcers.

- ◆ **Accurate:** 89.10% sensitivity, 95.90% specificity compared with PCR
- ◆ **Fast:** only 15 minutes required
- ◆ **Easy-to-use:** no special instruments required
- ◆ **Room temperature storage**

Cat.	Product
500070	StrongStep® HSV 1/2 Antigen Rapid Test

Detection	Specimen
Antigen	Mucocutaneous lesions swab