

StrongStep® Trichomonas/Candida Antigen Combo Rapid Test

An immunochromatographic assay for qualitative presumptive detection of *Trichomonas vaginalis* and/or *Candida albicans* antigens from vaginal swabs or wet mount saline solution



One time sampling and testing to get two results!

Background: *Trichomonas* infection is one of the most common, non-viral sexually transmitted disease worldwide, which is caused by a tiny single-celled parasite called *Trichomonas vaginalis*. Also, vulvovaginal candidiasis (VVC) is thought to be one of the most common causes of vaginal symptoms. About 20% of non-pregnant women aged 15 to 55 harbour *Candida albicans* in the vagina.

Vaginitis is one of the principal motives that lead women to seek out an obstetrician or gynecologist. Candidiasis and trichomoniasis are responsible for most of the cases of infectious vaginitis. Both of them need different treatment.

How can you benefit from StrongStep® Trichomonas/Candida Antigen combo rapid test

- ◆ **Fast:** Only 15 minutes required
- ◆ **Superior performance:** 93.60% sensitivity; 99.17% specificity for *Trichomonas* rapid test and 87.33% sensitivity; 99.26% specificity for *Candida albicans* rapid test compared with wet mount and culture
- ◆ **Save your time and cost:** One test for two pathogens with single swab
- ◆ **User-friendly:** Easily performed and interpreted by all health care personnel or even by patients themselves
- ◆ **Room temperature storage**

Cat.	Product
500040	<i>Trichomonas vaginalis</i>
500030	<i>Candida albicans</i>
500060	<i>Trichomonas/Candida</i> Combo

Detection	Specimen
Antigen	Vaginal discharge
Antigen	Vaginal discharge
Antigen	Vaginal discharge