



LIMING BIO

diagnostics are ASSURED

www.stddiagnostics.com www.limingbio.com

*Specialized In Manufacturing
Rapid Tests For
Sexually Transmitted Diseases*

StrongStep® Trichomonas/Candida Antigen Combo Rapid Test



An immunochromatographic assay for qualitative presumptive detection of Trichomonas vaginalis/Candida albicans antigen from vaginal discharge

- ◆ **Fast:** only 15 minutes required
- ◆ **Superior performance:** 93.60% sensitivity, 99.17% specificity for Trichomonas, 87.33% sensitivity, 99.26% specificity for Candida compared with wet mount and culture
- ◆ **Samples:** vaginal discharge, can be self-collected by patients
- ◆ **User-friendly:** easily performed and interpreted by all health care personnel
- ◆ **Room temperature storage**

Cat.	Product
500040	StrongStep® Trichomonas vaginalis Antigen Rapid Test
500030	Strong Step® Candida albicans Antigen Rapid Test
500060	Strong Step® Trichomonas/Candida Antigen Combo Rapid Test

Detection	Specimen
Antigen	Vaginal discharge
Antigen	Vaginal discharge
Antigen	Vaginal discharge